

TailorCut Gel Extraction Tool Set

Features


- Includes Gel Cutter, Gel Breaker and Gel Filter
- Easy and precise band excision
- Patent pending designed Gel Breaker fragments gel
- Increase Efficiency of Gel Purification
 - Enable precise excision of target band
 - Facilitate sample extraction from gel mix

Gel Cutter Tip



- Cutting size: 2 x 6mm
- Compatible with 1mL micropipette
- DNase and RNase free

Gel Breaker Tube




- Patent pending design fragments agarose and acrylamide gel slices
- Fragmented gels could be transferred using 1 mL pipet tips
- DNase and RNase free

Gel Filter Tube



- 5µm filter cartridge with 2mL microcentrifuge collection tube
- Maximum Force: 10,000xG
- DNase and RNase free



Order Online:

www.seqmatic.com
E-Mail: info@seqmatic.com
Tel: 510.870.0965

Catalog #	Product
TC-025	TailorCut Gel Extraction Set (25 Sets)